

R ONLY

DEVICE DESCRIPTION

The iFuse Implant System consists of cannulated triangular, titanium (Ti 6Al4V ELI) implants with a porous surface and an instrument system. Implant surface and shape are designed to prevent rotation and motion of the sacroiliac (SI) joint including acute, non-acute, and non-traumatic fractures involving the SI joint. The instrument system uses guide pins for accurate placement.

INDICATIONS

The iFuse Implant System is intended for sacroiliac fusion for the following conditions:

- Sacroiliac joint dysfunction that is a direct result of sacroiliac joint disruption and degenerative sacroiliitis. This includes conditions whose symptoms began during pregnancy or in the peripartum period and have persisted postpartum for more than 6 months.
- To augment immobilization and stabilization of the sacroiliac joint in skeletally mature patients undergoing sacropelvic fixation as part of a lumbar or thoracolumbar fusion.
- Acute, non-acute, and non-traumatic fractures involving the sacroiliac joint.

CONTRAINDICATIONS

1. Deformities or anatomic variations that prevent or interfere with iFuse placement.
2. Tumor of sacral or ilial bone.
3. Active infection at treatment site.
4. Unstable fracture of sacrum and/or ilium involving the sacroiliac joint.
5. Allergy to metal components.

WARNINGS

1. Women of childbearing potential should be cautioned that vaginal delivery of a fetus may not be advisable following SI joint fusion. If pregnancy occurs, the woman should review delivery options with her obstetrician.

PRECAUTIONS

1. Carefully read and follow all instructions prior to use.
2. Patient adherence to post-operative physical activity instructions is important to support long-term service life of the implant.
3. Pay careful attention to selection of implant size. Pre-operative X-rays and/or CT scan may be helpful in selecting optimal implant size.
4. Appropriate patient selection is necessary as patient factors such as size and weight may make use of the iFuse Implant System more difficult or impossible.
5. Inspect implants and instruments for damage prior to use. Do not use if damaged or worn. Do not attempt to repair.
6. Do not use any component from an opened or damaged

package.

7. Do not use implants after the expiration date.
8. If placing the implants in conjunction with an open procedure, the surgeon should take care not to destabilize the joint prior to placing the implants.

MRI SAFETY INFORMATION



Non-clinical testing and MRI simulations were performed to evaluate the entire family of the iFuse implant. Non-clinical testing demonstrated that the entire family of the iFuse implants is MR Conditional. A patient with implant(s) from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the implant(s) are expected to produce a maximum temperature rise of less than 5.4 °C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant(s) extends approximately 20 mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

RISKS

As with other surgical procedures used to treat SI joint conditions, including acute, non-acute, and non-traumatic fractures involving the SI joint, the risks associated with the iFuse Implant System surgical procedure include, but are not limited to the following:

1. Adverse reactions to anesthesia
2. Hemorrhage
3. Muscle damage
4. Hematoma or seroma
5. Neurological deficit, nerve root or peripheral nerve injury, irritation or damage
6. Vascular injury or damage that may result in catastrophic or fatal bleeding
7. Neurovascular injury
8. Damage to lymphatic vessels and/or lymphatic fluid exudation
9. Injury to intra-pelvic structures
10. Infection of the wound, deep infection, peritonitis
11. Wound dehiscence
12. Pulmonary or systemic embolism
13. Thrombosis, thrombophlebitis
14. Death
15. Bruising
16. Local swelling

17. Radiation exposure
18. Loss of fracture fixation

Potential risks specifically associated with the iFuse Implant System include, but are not limited to the following:

1. Infection
2. Pain, discomfort, or abnormal sensations due to presence of the implant
3. Instrument failure resulting in a complication
4. Migration, loosening or fracture of the implant
5. Pain in muscle(s) due to altered biomechanics
6. Nerve root or peripheral nerve irritation due to local swelling or altered biomechanics
7. Loss of fixation / stabilization
8. Metal sensitivity or allergic reaction
9. Failure to improve symptoms and/or function
10. Increased pain at treated or adjacent levels
11. Need for re-operation or removal of the implant(s)
12. Implant rejection
13. Response to wear debris
14. Decrease in bone density due to stress shielding
15. Failure to achieve SI joint fusion
16. Potential difficulty in delivering fetus vaginally due to device-related restriction of SI joint motion
17. Fracture of the sacrum and/or ilium

HOW SUPPLIED

Implants are provided sterile; do not resterilize. The instrumentation is provided separately, non-sterile, and must be sterilized prior to use following the *SI-BONE iFuse Family of Instruments Hospital Cleaning and Sterilization Instructions, U.S.A (P/N 300065)*.

STORAGE/HANDLING

1. Store packaged implants at room temperature.
2. Handle the iFuse Implant with care to prevent damage to the surface finish.

DIRECTIONS FOR USE

1. A minimum of two (2) implants per sacroiliac joint are recommended for treatment. Treatment typically involves at least three (3) iFuse Implants per SI joint.
2. For detailed information, refer to the relevant Surgical Technique Manual(s) prior to use of the iFuse Implant System.

GRAPHIC SYMBOLS GLOSSARY

	Caution, refer to instructions for use
	Sterilized using irradiation
	Non-sterile, sterilize prior to use
	Single sterile barrier system with protective packaging inside
	Batch code
	Use-by date
	Do not re-use
	Catalog number
	Do not use if package is damaged and consult instructions for use
	Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use or electronic instructions for use & symbols glossary www.si-bone.com/label
	MR Conditional
	Manufacturer
	Date of manufacture
	Medical device

 Manufactured for:
SI-BONE, Inc.
471 El Camino Real,
Suite 101
Santa Clara, CA 95050 USA

Customer Service:
USA: 408-207-0700 or Toll Free: 855-884-3873

Patents: www.si-bone.com