

Endovascular Today - US Device Guide - Carotid Artery Stenting Systems

Company Name	Product Name	Tapered Stent Diameters (Proximal/Distal) (mm)	Tapered Stent Lengths (mm)	Straight Stent Diameters (mm)	Straight Stent Lengths (mm)	Embolic Protection Device	Type	Position	US FDA Indicated Use
Abbott	RX Acculink Carotid Stent System	8/6, 10/7	30, 40	5, 6, 7, 8, 9, 10	20, 30, 40	Emboshield Nav6	BareWire filter delivery wire	Distal	High surgical risk: used in conjunction with Abbott Vascular’s Emboshield family of EPS, is indicated for the treatment of patients at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet the criteria outlined as follows: (1) patients with neurological symptoms and ≥ 50% stenosis of the common or internal carotid artery by ultrasound or angiogram or patients without neurological symptoms and ≥ 80% stenosis of the common or internal carotid artery by ultrasound or angiogram; (2) patients must have a reference vessel diameter within the range of 4–9 mm at the target lesion; standard surgical risk: used in conjunction with Abbott Vascular's Emboshield family of EPS, is indicated for the treatment of patients at standard risk for adverse events from carotid endarterectomy who require carotid revascularization and meet the criteria outlined as follows: (1) patients with neurological symptoms and ≥ 70% stenosis of the common or internal carotid artery by ultrasound or ≥ 50% stenosis of the common or internal carotid artery by angiogram or patients without neurological symptoms ≥ 70% stenosis of the common or internal carotid artery by ultrasound or ≥ 60% stenosis of the common or internal carotid artery by angiogram; (2) patients must have a reference vessel diameter within the range of 4–9 mm at the target lesion
Abbott	Xact Carotid Stent System	8/6, 9/7, 10/8	30, 40	7, 8, 9, 10	20, 30	Emboshield Nav6	BareWire filter delivery wire	Distal	The Xact Carotid Stent System, used in conjunction with the Emboshield family of Embolic Protection System is indicated for the improvement of the lumen diameter of carotid arteries in patients considered at high risk for adverse events from carotid endarterectomy who require percutaneous carotid angioplasty and stenting for occlusive artery disease and meet the criteria outlined below: patients with carotid artery stenosis (≥ 50% for symptomatic patients by ultrasound or angiography or ≥ 80% for asymptomatic patients by ultrasound or angiography), located between the origin of the common carotid artery and the intracranial segment of the internal carotid artery and patients must have a reference vessel diameter ranging between 4.8 mm and 9.1 mm at the target lesion

Company Name	Product Name	Tapered Stent Diameters (Proximal/Distal) (mm)	Tapered Stent Lengths (mm)	Straight Stent Diameters (mm)	Straight Stent Lengths (mm)	Embolic Protection Device	Type	Position	US FDA Indicated Use
Abbott	Xact™ Transcarotid Stent System	8/6, 9/7, 10/8	30, 40	7, 8, 9, 10	20, 30	Emboshield Nav6	BareWire filter delivery wire	Distal	The Xact Transcarotid Stent System, used in conjunction with the Enroute Transcarotid Neuroprotection System, is indicated for the improvement of the lumen diameter of carotid arteries in patients considered at high risk for adverse events from carotid endarterectomy who require percutaneous carotid angioplasty and stenting for occlusive artery disease and meet the criteria outlined below: patients with carotid artery stenosis ($\geq 50\%$ for symptomatic patients by ultrasound or angiography or $\geq 80\%$ for asymptomatic patients by ultrasound or angiography), located between the origin of the common carotid artery and the intracranial segment of the internal carotid artery and patients must have a reference vessel diameter ranging between 4.8 mm and 9.1 mm at the target lesion
Boston Scientific Corporation	Carotid Wallstent Endoprosthesis	–	–	6, 8, 10 (unconstrained)	21, 22, 24, 29, 31, 36, 37 (unconstrained)	FilterWire EZ Embolic Protection System	Filter	Distal	Used in conjunction with the Boston Scientific EPS, is indicated for the treatment of patients at high risk for adverse events from carotid endarterectomy due to either anatomic or comorbid conditions who require carotid revascularization in the treatment of ipsilateral or bilateral carotid artery disease and meet the criteria outlined: patients with neurological symptoms and $\geq 50\%$ stenosis of the common, internal carotid artery, and/or the bifurcation by ultrasound or angiography or patients without neurological symptoms and $\geq 80\%$ stenosis of the common, internal carotid artery, and/or the bifurcation by ultrasound or angiogram, and patients with a reference vessel diameter within the range of 4–9 mm at the target lesion
Boston Scientific Corporation	Enroute Transcarotid Stent System	8/6, 9/7, 10/8	30, 40	6, 7, 8, 9, 10	20 (6 mm only), 30, 40	Enroute Transcarotid Neuroprotection System	Robust flow reversal	Proximal	Used in conjunction with the Enroute Transcarotid Neuroprotection System (NPS), is indicated for the treatment of patients at high risk and standard risk for adverse events from carotid endarterectomy, who require carotid revascularization and meet the criteria outlined below: High-risk patients with neurological symptoms and $\geq 50\%$ stenosis of the common or internal carotid artery by ultrasound or angiogram or patients without neurological symptoms and $\geq 80\%$ stenosis of the common or internal carotid artery by ultrasound or angiogram; standard risk patients with neurological symptoms and $\geq 70\%$ stenosis of the common or internal carotid artery by ultrasound or $\geq 50\%$ stenosis of the common or internal carotid artery by angiogram or patients without neurological symptoms and $\geq$

Company Name	Product Name	Tapered Stent Diameters (Proximal/Distal) (mm)	Tapered Stent Lengths (mm)	Straight Stent Diameters (mm)	Straight Stent Lengths (mm)	Embolic Protection Device	Type	Position	US FDA Indicated Use
									70% stenosis of the common or internal carotid artery by ultrasound or ≥ 60% stenosis of the common or internal carotid artery by angiogram; all patients must have a vessel diameter of 4 to 9 mm at the target lesion, and carotid bifurcation must be located at minimum 5 cm above the clavicle to allow for placement of the Enroute Transcarotid NPS; please refer to instructions for use for indications, contraindications, warnings, and precautions
Contego Medical LLC	Neuroguard IEP System 3-in-1: stent + balloon + filter	9/7/8, 8/6/7 (hourglass shape)	30, 40	–	–	Integrated embolic protection (IEP) filter	Integrated embolic protection (IEP) filter	Distal	Indicated for improving the carotid luminal diameter in subjects at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet the following criteria: patients with symptomatic stenosis of the common or internal carotid artery with ≥ 50% as determined by angiography using NASCET methodology, or patients with asymptomatic stenosis of the common or internal carotid artery with ≥ 80% as determined by angiography using NASCET methodology; and patients with reference vessel diameters 4–8 mm; also indicated for postdilation of the stent component with simultaneous capture and removal of embolic material; the Neuroguard IEP System should always be used in conjunction with an available primary distal embolic protection device as described in the IFU
Cordis	Precise Nitinol Stent System	Autotapering	20, 30, 40	5, 6, 7, 8, 9, 10	20, 30, 40	Angioguard XP Emboli Capture Guidewire	Filter	Distal	Used in conjunction with the Cordis Angioguard XP Emboli Capture Guidewire, is indicated for treatment of patients at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet the following criteria: patients with neurological symptoms and ≥ 50% stenosis of the common or internal carotid artery by ultrasound or angiography or patients without neurological symptoms and ≥ 80% stenosis of the common or internal carotid artery by ultrasound or angiography and patients must have a vessel diameter of 4–9 mm at the target lesion; the vessel distal to the target lesion must be within the range of 3–7.5 mm to allow for placement of Cordis Angioguard XP Emboli Capture Guidewire
Cordis	Precise Pro RX Carotid Stent System	Autotapering	20, 30, 40	5, 6, 7, 8, 9, 10	20, 30, 40	Angioguard RX Emboli Capture Guidewire	Filter	Distal	Used in conjunction with the Cordis Angioguard RX Emboli Capture Guidewire, is indicated for treatment of patients at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet the following criteria:

Company Name	Product Name	Tapered Stent Diameters (Proximal/Distal) (mm)	Tapered Stent Lengths (mm)	Straight Stent Diameters (mm)	Straight Stent Lengths (mm)	Embollic Protection Device	Type	Position	US FDA Indicated Use
									patients with neurological symptoms and $\geq 50\%$ stenosis of the common or internal carotid artery by ultrasound or angiography or patients without neurological symptoms and $\geq 80\%$ stenosis of the common or internal carotid artery by ultrasound or angiography and patients must have a vessel diameter of 4–9 mm at the target lesion; the vessel distal to the target lesion must be within the range of 3–7.5 mm to allow for placement of Cordis Angioguard RX Guidewire
InspireMD	CGuard Prime	SmartFit auto-tapering	30, 40	8, 9, 10	30, 40	Emboshield Nav 6 or Mo.Ma or both	BareWire filter delivery wire (Nav6) or proximal embolic protection device (Mo.Ma)	Distal (Nav6) or proximal (Mo.Ma) or both	Used in conjunction with embolic protection devices specified in the labeling, is indicated for improving carotid luminal diameter in patients at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet both criteria outlined below: patients with neurological symptoms and $\geq 50\%$ stenosis of the common or internal carotid artery by angiogram, or patients without neurological symptoms and $\geq 80\%$ stenosis of the common or internal carotid artery by angiogram; patients having a vessel with reference diameters between 6.4 mm and 9 mm at the target lesion; the CGuard Prime Carotid Stent System should be used in conjunction with the Abbott Emboshield NAV6 Embolic Protection System or the Medtronic Mo.Ma Ultra Proximal Cerebral Protection Device
Medtronic	Protégé RX Carotid Stent System	8/6, 10/7	30, 40	6, 7, 8, 9, 10	20, 30, 40, 60	SpiderFX Embolic Protection Device	Filter: guidewire of choice	Distal	Used in conjunction with the SpiderFX embolic protection device, is indicated for the treatment of patients at high risk for adverse events from carotid endarterectomy who require percutaneous carotid revascularization and meet the following criteria: patients with carotid artery stenosis ($\geq 50\%$ for symptomatic patients by ultrasound or angiography or $\geq 80\%$ for asymptomatic patients by ultrasound or angiography) of the common or internal carotid artery, and patients must have a reference vessel diameter within the range of 4.5–9.5 mm at the target lesion
Terumo Interventional Systems	Roadsaver Carotid Artery Stent System	Self-tapering	Self-tapering	5, 6, 7, 8, 9, 10 (unconstrained)	22, 25, 33, 35, 37, 40, 47 (unconstrained)	Nanoparasol embolic protection system	Filter	Distal	When used in conjunction with the Nanoparasol Embolic Protection System (EPS), Roadsaver is indicated for the treatment of carotid artery stenosis in patients with elevated risk for adverse events following carotid endarterectomy and who meet the following criteria: 1) Patients who have either de novo atherosclerotic or post endarterectomy restenotic lesion(s) in

Company Name	Product Name	Tapered Stent Diameters (Proximal/Distal) (mm)	Tapered Stent Lengths (mm)	Straight Stent Diameters (mm)	Straight Stent Lengths (mm)	Embolic Protection Device	Type	Position	US FDA Indicated Use
									the internal carotid arteries or at the carotid bifurcation with ≥ 50% stenosis if symptomatic or ≥ 80% stenosis if asymptomatic (both defined by angiography), AND 2) patients having a vessel with reference diameters between 3.5–9.0 mm at the target lesion