

Endovascular Today - US Device Guide - Self-Expanding Stents								
Company Name	Product Name	Material Used	Introducer Size (F)	Guidewire Compatibility (inch)	Stent Diameter (mm)	Stent Length (mm)	Delivery System Length (cm)	US FDA Indicated Use
Abbott	Absolute Pro Vascular Stent	Nitinol	6	0.035	6, 7, 8, 9, 10	20, 30, 40, 60, 80, 100	80, 135	Iliac
Abbott	Supera Peripheral Stent System	Nitinol	6	0.018	4.5, 5, 5.5, 6, 6.5, 7.0, 7.5	20, 30, 40, 60, 80, 100, 120, 150	120	SFA/proximal popliteal
BD Interventional	E-Luminexx Biliary Stent	Nitinol	6	0.035	4, 5, 6, 7, 8, 9, 10, 12, 14	20, 30, 40, 60, 80, 100, 120	80, 135	Biliary
BD Interventional	E-Luminexx Vascular Stent	Nitinol	6	0.035	7, 8, 9, 10	20, 30, 40, 60, 80, 100	80, 135	Iliac
BD Interventional	LifeStar Biliary Stent System	Nitinol	6	0.035	6, 7, 8, 9, 10, 12, 14	20, 30, 40, 60, 80, 100, 120	80, 135	Biliary
BD Interventional	LifeStar Vascular Stent System	Nitinol	6	0.035	7, 8, 9, 10	20, 30, 40, 60, 80, 100	80, 135	Iliac
BD Interventional	LifeStent 5F Vascular Stent System	Nitinol	5	0.014–0.035	5, 6, 7	20, 30, 40, 60, 80, 100, 120, 150*, 170† (*only available in 5–6 mm diameters, †only available in 5 mm diameter)	80, 135	SFA, full popliteal
BD Interventional	LifeStent and LifeStent XL Biliary Stent System	Nitinol	6	0.035	5, 6, 7, 8, 9, 10	20, 30, 40, 60, 80, 100,* 120,* 150† (*only available in 5–7 mm diameters; †only available in 6, 7 mm diameters)	80, 130	Biliary
BD Interventional	LifeStent and LifeStent XL Vascular Stent System	Nitinol	6	0.035	5, 6, 7	20, 30, 40, 60, 80, 100, 120, 150, 170	80, 130	SFA, full popliteal
BD Interventional	LifeStent Solo Vascular Stent System	Nitinol	6	0.035	6, 7	200	100, 135	SFA, full popliteal
Boston Scientific Corporation	Epic Vascular Self-Expanding Stent System	Nitinol	6	0.035	6, 7, 8, 9, 10, 12	20, 30, 40, 60, 80, 100, 120	75, 120	Iliac
Boston Scientific Corporation	Innova Vascular Self-Expanding Stent System	Nitinol	6	0.035	5, 6, 7, 8	20, 40, 60, 80, 100, 120, 150, 200	75, 130	SFA/proximal popliteal
Boston Scientific	Venous Wallstent	Elgiloy	7–9	0.035	10, 12, 14, 16, 18, 20	20, 40, 42, 55, 60, 68, 80, 90, 94	75	Iliofemoral veins

Company Name	Product Name	Material Used	Introducer Size (F)	Guidewire Compatibility (inch)	Stent Diameter (mm)	Stent Length (mm)	Delivery System Length (cm)	US FDA Indicated Use
Corporation								
Boston Scientific Corporation	WallFlex Biliary Transhepatic Stent System	Platinol wire braid, Permalume covering	8, 8.5	0.035	8, 10	40, 60, 80, 100	75	Biliary
Boston Scientific Corporation	Wallstent Endoprosthesis	Elgiloy	6–10	0.035	10, 12, 14, 16	20, 40, 42, 60, 68, 90, 94	75	Central venous
Boston Scientific Corporation	Wallstent Endoprosthesis	Elgiloy	7–9	0.035	10, 12	40, 42, 60, 68, 90, 94	75	TIPS
Boston Scientific Corporation	Wallstent Endoprosthesis	Elgiloy	6–9	0.035	8, 10, 12	20, 40, 42, 60, 68, 80, 90, 94	75	Transhepatic biliary
Boston Scientific Corporation	Wallstent Endoprosthesis	Elgiloy	6–12	0.035	5, 6, 7, 8, 10, 12, 14, 16, 18, 20, 22, 24	20, 35, 40, 42, 45, 55, 60, 68, 70, 80, 90, 94	75, 135 (up to 12-mm diameter), 75 (14-mm diameter)	Tracheobronchial
Boston Scientific Corporation	Wallstent Endoprosthesis	Elgiloy	6	0.035	6, 7, 8, 9, 10	18, 20, 23, 24, 34, 35, 36, 38, 39, 46, 47, 49, 52, 55, 59, 61, 66, 67, 69	75, 135	Iliac
Cook Medical	Zilver 518 Biliary Stent	Nitinol	5	0.018	4, 5	20, 30, 40, 60, 80	125	Biliary
Cook Medical	Zilver 518 Vascular Stent	Nitinol	5	0.018	6, 7, 8, 9, 10	20, 30, 40, 60, 80	125	Iliac
Cook Medical	Zilver 635 Biliary Stent	Nitinol	6	0.035	4, 5, 12, 14	4-mm diameter: 80-mm length; 5-mm diameter: 20-, 30-, 40-, 60-, 80-mm lengths; 8-mm diameter: 40-, 60-, 80-mm lengths; 10-mm diameter: 40-, 60-, 80-mm lengths; 12-, 14-mm diameter: 30-, 40-, 60-, 80-mm lengths	40, 80, 125	Biliary
Cook Medical	Zilver 635 Vascular Stent	Nitinol	6	0.035	6, 7, 8, 9, 10	20, 30, 40, 60, 80	80, 125	Iliac
Cook Medical	Zilver Flex 35 Biliary Stent	Nitinol	6	0.035	5, 6, 7, 8, 9, 10	100, 120, 140	80, 125	Biliary

Company Name	Product Name	Material Used	Introducer Size (F)	Guidewire Compatibility (inch)	Stent Diameter (mm)	Stent Length (mm)	Delivery System Length (cm)	US FDA Indicated Use
Cordis	Smart Control Biliary Stent System	Nitinol	6	0.035	12, 14	30, 40, 60, 80	80, 120	Biliary
Cordis	Smart Control Iliac Stent System	Nitinol	6	0.035	9, 10	20, 30, 40, 60	80, 120	Iliac
Cordis	Smart Control Vascular Stent System	Nitinol	6	0.035	6, 7, 8	20, 30, 40, 60, 80, 100	80, 120	Iliac, SFA
Cordis	Smart Flex Biliary Stent System	Nitinol	6	0.035	5, 6, 7, 8, 9, 10	30, 40, 60, 80, 100, 120, 150	80, 120	Biliary
Cordis	Smart Vascular Stent System	Nitinol	6	0.035	6, 7, 8	120, 150	120	SFA
Medtronic	EverFlex Self-Expanding Peripheral Stent System	Nitinol	6	0.035	6, 7, 8	20, 30, 40, 60, 80, 100, 120	80, 120	Iliac
Medtronic	EverFlex Self-Expanding Peripheral Stent System	Nitinol	6	0.035	6, 7, 8	20, 30, 40, 60, 80, 100, 120, 150, 200	80, 120	SFA, PPA
Medtronic	EverFlex Self-Expanding Peripheral Stent With Entrust Delivery System	Nitinol	5	0.035	6, 7, 8	20, 40, 60, 80, 100, 120, 150	80, 120, 150	SFA/PPA
Medtronic	Protege EverFlex Self-Expanding Stent System	Nitinol	6	0.035	5	20, 30, 40, 60, 80, 100, 120	80, 120	Biliary
Medtronic	Protege GPS Self-Expanding Peripheral Stent System	Nitinol	6	0.035	9, 10, 12	20, 30, 40, 60, 80	80, 120	Iliac
Medtronic	Protege GPS Self-Expanding Stent System	Nitinol	6	0.035	9, 10, 12, 14	20, 30, 40, 60, 80	80, 120	Biliary
Teleflex	Astron	Nitinol	6	0.035	7, 8, 9, 10	30, 40, 60, 80	72, 130	Indicated for improving luminal diameter in patients with iliac atherosclerotic lesions in vessel reference diameters between 4.3–9.5 mm and lesion lengths up to 105 mm
Teleflex	Pulsar-18	Nitinol	4	0.018	4–7	20–150	90, 135	Indicated for use to improve luminal diameter in patients with symptomatic de novo, restenotic, or occlusive lesions located in the

Company Name	Product Name	Material Used	Introducer Size (F)	Guidewire Compatibility (inch)	Stent Diameter (mm)	Stent Length (mm)	Delivery System Length (cm)	US FDA Indicated Use
								superficial femoral or proximal popliteal arteries, with reference vessel diameters from 3–6 mm and total lesion lengths up to 190 mm
Teleflex	Pulsar-18 T3	Nitinol	4	0.018	4, 5, 6, 7	20, 30, 40, 60, 80, 100, 120, 150, 170, 200	90, 135	Indicated for use to improve luminal diameter in patients with symptomatic de novo, restenotic or occlusive lesions located in the superficial femoral or proximal popliteal arteries, with reference vessel diameters from 3 to 6 mm and total lesion lengths up to 190 mm
Terumo Interventional Systems	R2P Misago RX Self-Expanding Peripheral Stent	Nitinol	6	0.035	6, 7, 8	40, 60, 80, 100, 120, 150	200	SFA, proximal popliteal
Veryan Medical	BioMimics 3D Vascular Stent System	Nitinol	6	0.035	5, 6, 7	60, 80, 100, 125, 150	113	SFA, proximal popliteal
Vesalio Inc.	NeVa VS - Vasospasm Treatment Device	Nitinol	Can be delivered through catheters with minimum ID of 0.021 inch (1.6 F)	Max: 0.018	4	Full length: 35 mm; working length: 30 mm (excluding the 5-mm tip)	200	Endovascular treatment of cerebral vasospasm after aneurysmal SAH